

Chromatography Solutions

Avantor® ACE® Generix®

High performance, excellent value
HPLC and UHPLC Columns



Avantor® ACE® Generix®

HIGH PERFORMANCE, EXCELLENT VALUE HPLC AND UHPLC COLUMNS

High Efficiency 1.8, 3 and 5 µm Particles

- High purity spherical silica
- Ultra-narrow particle size distribution for excellent column efficiency
- C18, C8, Phenyl-Hexyl and SIL phases

Excellent Reproducibility

- Manufactured in the UK in ISO9001 and ISO14001 accredited facilities
- Sharp, highly symmetrical peaks for acidic, basic and neutral analytes
- Every column individually tested and supplied with a test chromatogram
- Excellent column-to-column and batch-to-batch reproducibility

Excellent Column Lifetime

- Highly stable bonded phases for long column lifetimes
- Ultra-robust HPLC and UHPLC column packing technology
- High performance guard system available



Phase Specifications

Bonded Phase ^a	Particle Size (µm)	Surface Area (m ² /g)	Pore Size (Å)	Carbon Load (%)	Recommended pH Range ^b	USP Listing
C18(2)	1.8, 3 and 5	350	95	20.3	2.0-8.0	L1
C8(2)	1.8, 3 and 5	350	95	11.9	2.0-8.0	L7
Phenyl-Hexyl	1.8, 3 and 5	350	95	15.2	2.0-8.0	L11
SIL	3 and 5	350	95	-	2.0-8.0	L3


^aThe full range of Avantor® ACE® HPLC and UHPLC columns is available with 1.7, 2, 3, 5 and 10 µm particle sizes and a wider range of stationary phase options.

^bFor optimum column lifetime, a pH range of 2-8 is recommended. To increase column lifetime at mid/high pH, organic buffers, low buffer concentrations, high % organic solvent and low temperatures must be considered. Further information is contained within "A Guide to HPLC and LC/MS Buffer Selection" by John Dolan.

SEPARATION OF VITAMINS

AVANTOR® ACE® GENERIX® 1.8 PHENYL-HEXYL


Column: Avantor® ACE® Generix® 1.8 Phenyl-Hexyl, 100 x 2.1 mm
Part Number: GEN-18PHEX-1002U
Mobile Phase: A: 20 mM KH₂PO₄ (pH 2.7) in H₂O
 B: 20 mM KH₂PO₄ (pH 2.7) in MeOH/H₂O, 50:50 v/v
Gradient: 5 to 95 % B in 10 minutes
Flow Rate: 0.2 mL/min
Temperature: 30°C
Detection: 210 nm
Sample: 1: Pyridoxine (B6); 2: Pantothenic acid (B5);
 3: p-Aminobenzoic acid; 4: Cyanocobalamin (B12);
 5: Biotin (B7)



SEPARATION OF FATTY ACIDS

AVANTOR® ACE® GENERIX® 5 C8(2)

Column: Avantor® ACE® Generix® 5 C8(2), 250 x 4.6 mm
Part Number: GEN-5C82-2546
Mobile Phase: MeCN/H₂O 90:10 v/v
Flow Rate: 1.0 mL/min
Temperature: 22°C
Detection: ELSD
Sample: 1: Linolenic acid; 2: Myristic acid; 3: Linoleic acid;
 4: Palmitic acid; 5: Oleic acid; 6: Stearic acid



Avantor® ACE® Generix®


- Avantor® ACE® Generix® has been specifically developed as a high quality, cost effective alternative for many leading and legacy HPLC column brands.
- Combining efficiency, reproducibility and good column lifetimes, Avantor® ACE® Generix® is an excellent option to help reduce your analysis costs.

SEPARATION OF PHARMACEUTICALLY RELEVANT COMPOUNDS


AVANTOR® ACE® GENERIX® CAN BE USED AS A COST EFFECTIVE ALTERNATIVE TO MANY LEADING AND LEGACY COLUMN BRANDS:

Column: 150 x 4.6 mm
Mobile Phase: A: 20 mM KH₂PO₄ (pH 2.7) in H₂O
 B: 20 mM KH₂PO₄ (pH 2.7) in 65:35 MeCN/H₂O v/v
Gradient: 5 to 100 % B in 18 minutes
Flow Rate: 1.0 mL/min
Temperature: 40°C
Detection: 214 nm
Sample: 1: Procainamide; 2: Amiloride; 3: Caffeine; 4: Methyl Phenyl Sulfoxide;
 5: Propiophenone


Avantor® ACE® Generix® 3 C18(2)




Phenomenex Luna 3 C18




Waters SunFire 3.5 C18




Waters Symmetry 3.5 C18



Agilent Zorbax 3.5 SB-C18



Agilent Zorbax Eclipse 3.5 XDB-C18



Comparative data above may not be representative of different applications. Analysts must verify suitability for specific requirements/analytes. Avantor®, ACE® and Generix® are registered trademarks of Avantor, Inc. or its affiliates. All other trademarks are the property of their respective owners.

Avantor® ACE® Generix® C18(2) Cross Reference Guide

For more information
or advice please
contact our technical
support team at
[chromsupport@
avantorsciences.com](mailto:chromsupport@avantorsciences.com)

ALTERNATIVES TO THE FOLLOWING LEADING AND LEGACY COLUMN BRANDS:

(Manufacturers' original part numbers are shown in the table below. Avantor® ACE® Generix® columns are also available in additional dimensions - please enquire for further details.)

Column Brand	Manufacturer	1.8 µm		3 µm		
		2.1 x 50 mm	2.1 x 100 mm	2.1 x 50 mm	4.6 x 100 mm	4.6 x 150 mm
Avantor® ACE® Generix® C18(2)		GEN-18C182-0502	GEN-18C182-1002	GEN-3C182-0502	GEN-3C182-1046	GEN-3C182-1546
Acclaim 120 C18	Thermo Scientific	-	-	059128	059132	059133
Ascentis C18	Supelco	-	-	581300-U	581321-U	581322-U
Betasil C18	Thermo Scientific	-	-	-	-	-
Discovery C18	Supelco	-	-	-	-	-
Discovery HS C18	Supelco	-	-	569253-U	568521-U	569252-U
Force C18	Restek	9634252	9634212	9634352	9634315	9634365
Fortis C18 ¹	Fortis Technologies	F18-0203-01	F18-0205-01	F18-0203-03	F18-0505-03	F18-0507-03
GraceSmart C18	Grace	-	-	5141756	5141753	5141752
Hypersil ODS	Thermo Scientific	-	-	-	30103-104630	30103-154630
Hypersil BDS C18	Thermo Scientific	-	-	28103-052130	28103-104630	28103-154630
Hypersil GOLD C18 Selectivity ²	Thermo Scientific	-	-	25003-052130	25003-104630	25003-154630
Hypersil HyPURITY C18	Thermo Scientific	-	-	22103-052130	22103-104630	22103-154630
Inertsil ODS	GL Sciences	-	-	-	-	-
Inertsil ODS-2	GL Sciences	-	-	-	-	-
Inertsil ODS-3 ³	GL Sciences	5020-84652	5020-84654	5020-04412	5020-01775	5020-01771
Inertsil ODS-4 ³	GL Sciences	5020-81202	5020-81204	5020-04012	5020-04044	5020-04045
InertSustain C18 ³	GL Sciences	5020-14352	5020-14354	5020-07412	5020-07444	5020-07445
LiChrosorb RP-18	Merck	-	-	-	-	-
LiChrospher RP-18	Merck	-	-	-	-	-
LiChrospher RP-18e	Merck	-	-	-	-	-
Luna C18(2) ¹	Phenomenex	-	-	-	-	00F-4251-E0
Nucleosil 100 C18	Macherey-Nagel	-	-	-	720841.46	720949.46
Nucleosil 120 C18	Macherey-Nagel	-	-	-	720149.46	720740.46
Nucleosil C18 HD	Macherey-Nagel	-	-	-	-	720193.46
Nucleodur C18 Gravity ⁵	Macherey-Nagel	760079.20	760076.20	760080.20	760084.46	760083.46
Nucleodur C18 ec ⁵	Macherey-Nagel	-	-	760050.20	760054.46	760053.46
Prodigy ODS-2	Phenomenex	-	-	-	-	-
Prodigy ODS-3 ⁵	Phenomenex	-	-	-	00D-4222-E0	00F-4222-E0
Prodigy ODS-3V	Phenomenex	-	-	-	-	00F-4243-E0
Purospher STAR RP-18e ³	Merck	1.50646.0001	1.50648.0001	1.50651.0001	1.51448.7184	1.50470.0001
Roc C18	Restek	-	-	-	9534315	9534365
SunFire C18 ⁴	Waters	-	-	186002533	186002553	186002554
Supelcosil LC-18	Supelco	-	-	-	-	58985
Supelcosil LC-18-DB	Supelco	-	-	-	-	58993
Synchronis C18 ¹	Thermo Scientific	97102-052130	97102-102130	-	97103-104630	97103-154630
Symmetry C18 ⁴	Waters	-	-	WAT200650	WAT066220	WAT200632
Ultra C18	Restek	-	-	9174352	9174315	9174365
Vision HT C18 HighLoad ^{6,6}	Grace	5142546	5142547	5151981	5151972	5151971
Waters Spherisorb ODS1	Waters	-	-	-	PSS833412	PSS833413
Waters Spherisorb ODS2	Waters	-	-	-	PSS832112	PSS832113
YMC-Pack ODS-A	YMC	-	-	AA12S03-05Q1QT	AA12S03-1046WT	AA12S03-1546WT
YMC-Pack ODS-AM	YMC	-	-	AM12S03-05Q1QT	AM12S03-1046WT	AM12S03-1546WT
YMC-Pack Pro C18 ^{3,5}	YMC	AS12S02-0502WT	-	AS12S03-05Q1QT	AS12S03-1046WT	AS12S03-1546WT
Zorbax Eclipse XDB-C18 ⁴	Agilent	927700-902	928700-902	971700-902	961967-902	963967-902
Zorbax Eclipse Plus C18 ⁴	Agilent	959757-902	959758-902	959743-902	959961-902	959963-902
Zorbax C18	Agilent	-	-	-	-	-
Zorbax Rx-C18 ⁴	Agilent	-	-	-	861967-902	863967-902
Zorbax SB-C18 ⁴	Agilent	857700-902	858700-902	871700-902	861953-902	863953-902

¹ 1.7 µm; ² 1.9 µm; ³ 2 µm; ⁴ 3.5 µm; ⁵ 2.0 mm ID; ⁶ 1.5 µm, ⁷ 3 µm

Avantor®, ACE® and Generix® are registered trademarks of Avantor, Inc. or its affiliates. All other trademarks are the property of their respective owners.

Column Brand	Manufacturer	5 µm				
		2.1 x 50 mm	4.6 x 50 mm	4.6 x 100 mm	4.6 x 150 mm	4.6 x 250 mm
Avantor® ACE® Generix® C18(2)		GEN-5C182-0502	GEN-5C182-0546	GEN-5C182-1046	GEN-5C182-1546	GEN-5C182-2546
Acclaim 120 C18	Thermo Scientific	059142	059146	059147	059148	059149
Ascentis C18	Supelco		581323-U	-	581324-U	581325-U
Betasil C18	Thermo Scientific	70105-052130	70105-054630	70105-1504630	70105-154630	70105-254630
Discovery C18	Supelco	50494721	504947	-	504955	504971
Discovery HS C18	Supelco	568500-U	568520-U	-	568522-U	568523-U
Force C18	Restek	9634552	-	9634515	9634565	9634575
Fortis C18 ¹	Fortis Technologies	F18-0203-05	F18-0503-05	F18-0505-05	F18-0507-05	F18-0509-05
GraceSmart C18	Grace	-	-	-	5138812	5138810
Hypersil ODS	Thermo Scientific	30105-052130	30105-054630	30105-104630	30105-154630	30105-254630
Hypersil BDS C18	Thermo Scientific	28105-052130	28105-054630	28105-104630	28105-154630	28105-254630
Hypersil GOLD C18 Selectivity ²	Thermo Scientific	25005-052130	25005-054630	25005-104630	25005-154630	25005-254630
Hypersil HyPURITY C18	Thermo Scientific	22105-052130	22105-054630	22105-104630	22105-154630	22105-254630
Inertsil ODS	GL Sciences	-	-	-	5020-02124	5020-02128
Inertsil ODS-2	GL Sciences	-	-	-	5020-01124	5020-01128
Inertsil ODS-3 ³	GL Sciences	5020-04512	5020-01763	5020-01765	5020-01731	5020-01732
Inertsil ODS-3V	GL Sciences	-	-	-	5020-01801	5020-01802
Inertsil ODS-4 ³	GL Sciences	5020-03912	5020-03942	5020-03944	5020-03945	5020-03946
Inertsil ODS-4V	GL Sciences	-	-	-	5020-10901	5020-10902
InertSustain C18 ³	GL Sciences	5020-07312	5020-07342	5020-07344	5020-07345	5020-07346
Luna C18(2) ⁵	Phenomenex	00B-4252-B0	00B-4252-E0	00D-4252-E0	00F-4252-E0	00G-4252-E0
LiChrosorb RP-18	Merck	-	-	-	1.50009.7052	1.51902.0001
LiChrospher RP-18	Merck	-	-	1.50600.0001	1.50601.0001	1.50602.0001
LiChrospher RP-18e	Merck	-	-	1.50603.0001	1.50604.0001	1.50605.0001
Nucleosil 100 C18	Macherey-Nagel	-	-	720141.46	720120.46	720014.46
Nucleosil 120 C18	Macherey-Nagel	-	-	-	720730.46	720041.46
Nucleosil C18 HD	Macherey-Nagel	-	-	-	720294.46	720280.46
Nucleodur C18 Gravity ⁵	Macherey-Nagel	760102.20	760102.46	760104.46	760103.46	760101.46
Nucleodur C18ec ⁵	Macherey-Nagel	760004.20	760004.46	760013.46	760008.46	760002.46
Prodigy ODS-2	Phenomenex	-	00B-3300-E0	00D-3300-E0	00F-3300-E0	00G-3300-E0
Prodigy ODS-3 ³	Phenomenex	00B-4097-B0	00B-4097-E0	00D-4097-E0	00F-4097-E0	00G-4097-E0
Prodigy ODS-3V	Phenomenex	-	-	-	00F-4241-E0	00G-4241-E0
Purospher STAR RP-18e	Merck	-	-	1.50622.0001	1.51455.0001	1.51456.0001
Roc C18	Restek	-	-	9534515	9534565	9534575
SunFire C18 ⁴	Waters	186002539	186002557	186002558	186002559	186002560
Supelcosil LC-18	Supelco	-	-	59209	58230-U	58298
Supelcosil LC-18-DB	Supelco	-	58345	59208	58348	58355-U
Synchronis C18 ¹	Thermo Scientific	97105-052130	97105-054630	97105-104630	97105-154630	97105-254630
Symmetry C18 ⁴	Waters	186000206	186000207	186002616	WAT045905	WAT054275
Ultra C18	Restek	9174552	9174555	9174515	9174565	9174575
Vision HT C18 HighLoad ⁶	Grace	-	-	-	5151970	5151920
Waters Spherisorb ODS1	Waters	-	PSS830611	PSS830612	PSS830613	PSS830615
Waters Spherisorb ODS2	Waters	-	PSS831911	PSS831912	PSS831913	PSS831915
YMC-Pack ODS-A	YMC	AA12S05-05Q1WT	AA12S05-0546WT	AA12S05-1046WT	AA12S05-1546WT	AA12S05-2546WT
YMC-Pack ODS-AM	YMC	AM12S05-0502WT	AM12S05-0546WT	AM12S05-1046WT	AM12S05-1546WT	AM12S05-2546WT
YMC-Pack Pro C18 ^{3,5}	YMC	AS12S05-0502WT	AS12S05-0546WT	AS12S05-1046WT	AS12S05-1546WT	AS12S05-2546WT
Zorbax Eclipse XDB-C18 ⁴	Agilent	960967-902	946975-902	-	993967-902	990967-902
Zorbax Eclipse Plus C18 ⁴	Agilent	959746-902	959946-902	959996-902	959993-902	959990-902
Zorbax ODS	Agilent	-	-	-	883952-702	880952-702
Zorbax Rx-C18 ⁴	Agilent	-	-	-	883967-902	880967-902
Zorbax SB-C18 ⁴	Agilent	860975-902	846975-902	-	883975-902	880975-902

¹ 1.7 µm; ² 1.9 µm; ³ 2 µm; ⁴ 3.5 µm; ⁵ 2.0 mm ID; ⁶ 1.5 µm

Avantor® ACE® Generix® C8(2) Cross Reference Guide

ALTERNATIVES TO THE FOLLOWING LEADING AND LEGACY COLUMN BRANDS:

(Manufacturers' original part numbers are shown in the table below.
Avantor® ACE® Generix® columns are also available in additional
dimensions - please enquire for further details.)

Column Brand	Manufacturer	1.8 µm		3 µm		
		2.1 x 50 mm	2.1 x 100 mm	2.1 x 50 mm	4.6 x 100 mm	4.6 x 150 mm
Avantor® ACE® Generix® C8(2)	Avantor			GEN-3C82-0502	GEN-3C82-1046	GEN-3C82-1546
Acclaim 120 C8	Thermo Scientific	-	-	059122	-	059127
Ascentis C8	Supelco	-	-	581400-U	581407-U	581408-U
Betasil C8	Thermo Scientific	-	-	-	-	-
Discovery C8	Supelco	-	-	-	-	-
Fortis C8 ¹	Fortis Technologies	F08-0203-01	F08-0205-01	F08-0203-03	F08-0505-03	F08-0507-03
Hypersil MOS	Thermo Scientific	-	-	-	-	-
Hypersil BDS C8	Thermo Scientific	-	-	28203-052130	28203-104630	28203-154630
Hypersil GOLD C8 ²	Thermo Scientific	25202-052130	25202-102130	25203-052130	25203-104630	25203-154630
Inertsil C8	GL Sciences	-	-	-	-	-
Inertsil C8-3 ³	GL Sciences	5020-84931	5020-84933	5020-04812	5020-04844	5020-01911
Inertsil C8-4 ³	GL Sciences	5020-81282	5020-81284	5020-03972	5020-03995	5020-03996
InertSustain C8 ³	GL Sciences	5020-16236	5020-16238	5020-16133	5020-16156	5020-16157
Luna C8(2) ⁴	Phenomenex	-	-	00B-4248-B0	00D-4248-E0	00F-4248-E0
LiChrosorb RP-8	Merck	-	-	-	-	-
LiChrospher RP-8	Merck	-	-	-	-	-
LiChrospher RP-8e	Merck	-	-	-	-	-
Nucleosil 100 C8	Macherey-Nagel	-	-	-	-	-
Nucleosil 120 C8	Macherey-Nagel	-	-	-	-	720214.46
Nucleosil C8 HD	Macherey-Nagel	-	-	-	-	-
Nucleodur C8 Gravity ⁴	Macherey-Nagel	760755.20	760757.20	760653.20	760656.46	760658.46
Nucleodur C8 ec ⁴	Macherey-Nagel	-	-	760063.20	760059.46	760061.46
Prodigy C8 ⁴	Phenomenex	-	-	-	-	-
Purospher STAR RP-8e ³	Merck	1.50630.0001	1.50629.0001	1.50674.0001	-	1.50009.7220
Roc C8	Restek	-	-	-	9533315	9533365
SunFire C8 ⁵	Waters	-	-	186002710	186002731	186002732
Supelcosil LC-8	Supelco	-	-	-	-	58983
Supelcosil LC-8 DB	Supelco	-	-	-	-	58991
Synchronis C8 ¹	Thermo Scientific	97202-052130	97202-102130	-	97203-104630	97203-154630
Symmetry C8 ⁵	Waters	-	-	WAT200624	WAT066204	WAT200630
Ultra C8	Restek	-	-	9103352	9103315	9103365
Waters Spherisorb C84	Waters	-	-	-	PSS832212	PSS832213
YMC-Pack C8	YMC	-	-	OC12S03-05Q1QT	OC12S03-1046WT	OC12S03-1546WT
YMC-Pack Pro C8	YMC	-	-	OS12S03-05Q1QT	OS12S03-1046WT	OS12S03-1546WT
Zorbax Eclipse XDB-C8 ⁵	Agilent	927700-906	928700-906	971700-906	961967-906	963967-906
Zorbax Eclipse Plus C8 ⁵	Agilent	959757-906	959758-906	959743-906	959961-906	959963-906
Zorbax Rx-C8 ⁵	Agilent	-	-	-	-	-
Zorbax SB-C8 ⁵	Agilent	857700-906	858700-906	871700-906	861953-906	863953-906

¹ 1.7 µm; ² 1.9 µm; ³ 2 µm; ⁴ 2.0 mm ID; ⁵ 3.5 µm

Avantor®, ACE® and Generix® are registered trademarks of Avantor, Inc. or its affiliates.
All other trademarks are the property of their respective owners



For more information or advice please contact our technical support team at chromsupport@avantorsciences.com

Column Brand	Manufacturer	5 µm				
		2.1 x 50 mm	4.6 x 50 mm	4.6 x 100 mm	4.6 x 150 mm	4.6 x 250 mm
Avantor® ACE® Generix® C8(2)	Avantor	GEN-5C82-0502	GEN-5C82-0546	GEN-5C82-1046	GEN-5C82-1546	GEN-5C82-2546
Acclaim 120 C8	Thermo Scientific	059134	059138	059139	059140	059141
Ascentis C8	Supelco	-	-	-	581424-U	581425-U
Betasil C8	Thermo Scientific	-	-	-	70205-154630	70205-254630
Discovery C8	Supelco	59352-U21	59352-U	569423-U	59353-U	59354-U
Fortis C8 ¹	Fortis Technologies	F08-0203-05	F08-0503-05	F08-0505-05	F08-0507-05	F08-0509-05
Hypersil MOS	Thermo Scientific	-	30205-054630	30205-104630	30205-154630	30205-254630
Hypersil BDS C8	Thermo Scientific	28205-052130	28205-054630	28205-104630	28205-154630	28205-254630
Hypersil GOLD C8 ²	Thermo Scientific	25205-052130	25205-054630	25205-104630	25205-154630	25205-254630
Inertsil C8 ³	GL Sciences	-	-	-	5020-01224	5020-01228
Inertsil C8-3 ³	GL Sciences	5020-04912	5020-04942	5020-04944	5020-01900	5020-01901
Inertsil C8-4 ³	GL Sciences	5020-04052	5020-04082	5020-04084	5020-04085	5020-04086
InertSustain C8	GL Sciences	5020-16003	5020-16024	5020-16026	5020-16027	5020-16028
Luna C8(2) ⁴	Phenomenex	00B-4249-B0	00B-4249-E0	00D-4249-E0	00F-4249-E0	00G-4249-E0
LiChrosorb RP-8	Merck	-	-	-	-	1.51903.0001
LiChrospher RP-8	Merck	-	-	1.50634.0001	1.50635.0001	1.50636.0001
LiChrospher RP-8e	Merck	-	-	1.50637.0001	1.50638.0001	1.50639.0001
Nucleosil 100 C8	Macherey-Nagel	-	-	-	720990.46	720103.46
Nucleosil 120 C8	Macherey-Nagel	-	-	-	720735.46	720052.46
Nucleosil C8 HD	Macherey-Nagel	-	-	-	720194.46	720196.46
Nucleodur C8 Gravity ⁴	Macherey-Nagel	760750.20	760750.46	760754.46	760752.46	760753.46
Nucleodur C8 ec ⁴	Macherey-Nagel	760700.20	760700.46	760704.46	760702.46	760703.46
Prodigy C8 ⁴	Phenomenex	-	00B-3301-E0	00D-3301-E0	00F-3301-E0	00G-3301-E0
Purospher STAR RP-8e	Merck	-	-	-	1.50031.0001	1.50032.0001
Roc C8	Restek	-	-	9533515	9533565	9533575
SunFire C8 ⁵	Waters	186002715	186002735	186002736	186002737	186002738
Supelcosil LC-8	Supelco	-	58238	-	58220-U	58297
Supelcosil LC-8-DB	Supelco	-	58344	-	58347	58354
Synchronis C8 ¹	Thermo Scientific	-	-	97205-104630	97205-154630	97205-254630
Symmetry C8 ⁵	Waters	186000212	186000213	186002617	WAT045995	WAT054270
Ultra C8	Restek	9103552	9103555	9103515	9103565	9103575
Waters Spherisorb C81	Waters	-	PSS831811	PSS831812	PSS831813	PSS831815
YMC-Pack C8	YMC	OC12S05-05Q1WT	OC12S05-0546WT	OC12S05-1046WT	OC12S05-1546WT	OC12S05-2546WT
YMC-Pack Pro C8 ⁴	YMC	OS12S05-0502WT	OS12S05-0546WT	OS12S05-1046WT	OS12S05-1546WT	OS12S05-2546WT
Zorbax Eclipse XDB-C8 ⁵	Agilent	960967-906	946975-906	-	993967-906	990967-906
Zorbax Eclipse Plus C8 ⁵	Agilent	959746-906	959946-906	959996-906	959993-906	959990-906
Zorbax Rx-C8	Agilent	-	-	-	883967-901	880967-901
Zorbax SB-C8 ⁵	Agilent	860975-906	846975-906	-	883975-906	880975-906

¹ 1.7 µm; ² 1.9 µm; ³ 2 µm; ⁴ 2.0 mm ID; ⁵ 3.5 µm

Avantor® ACE® Generix® Phenyl-Hexyl Cross Reference Guide

ALTERNATIVES TO THE FOLLOWING LEADING AND LEGACY COLUMN BRANDS:

(Manufacturers' original part numbers are shown in the table below. Avantor® ACE® Generix® columns are also available in additional dimensions - please enquire for further details.)

Column Brand	Manufacturer	1.8 µm		3 µm		
		2.1 x 50 mm	2.1 x 100 mm	2.1 x 50 mm	4.6 x 100 mm	4.6 x 150 mm
Avantor® ACE® Generix® Phenyl Hexyl		GEN-18PHEX-0502	GEN-18PHEX-1002	GEN-3PHEX-0502	GEN-3PHEX-1046	GEN-3PHEX-1546
Acclaim 120 Phenyl-1	Thermo Scientific	-	-	-	-	071969
Hypersil Phenyl	Thermo Scientific	-	-	-	-	-
Hypersil Phenyl-BES	Thermo Scientific	-	-	-	-	28903-154630
Hypersil GOLD Phenyl ¹	Thermo Scientific	25902-052130	25902-102130	-	25903-104630	25903-154630
Inertsil Ph	GL Sciences	-	-	-	-	-
Inertsil Ph-3 ²	GL Sciences	5020-85131	5020-85133	5020-05012	5020-05044	5020-01931
InertSustain Phenyl ²	GL Sciences	5020-16536	5020-16538	5020-16433	5020-16456	5020-16457
InertSustain Phenylhexyl	GL Sciences	-	-	5020-89125	5020-89148	5020-89149
Luna Phenyl-Hexyl ³	Phenomenex	-	-	00B-4256-B0	00D-4256-E0	00F-4256-E0
Nucleosil C6H5	Macherey-Nagel	-	-	-	-	-
Nucleodur Phenyl-Hexyl ³	Macherey-Nagel	760563.20	760566.20	760573.20	760576.46	760578.46
Prodigy Phenyl-3 ³	Phenomenex	-	-	-	-	-
Purospher STAR Phenyl ²	Merck	1.51013.0001	1.51014.0001	1.50672.0001	-	-
Roc Phenyl-Hexyl	Restek	-	-	-	9535315	9535365
Synchronis Phenyl ⁴	Thermo Scientific	97902-052130	97902-102130	-	-	-
Waters Spherisorb Phenyl	Waters	-	-	-	PSS833712	PSS833713
YMC-Pack Ph ³	YMC	-	-	PH12S03-0502WT	PH12S03-1046WT	PH12S03-1546WT
Zorbax Eclipse XDB-Phenyl ⁵	Agilent	-	-	-	-	963967-912
Zorbax Eclipse Plus-Phenyl-Hexyl ⁵	Agilent	959757-912	959758-912	959743-912	959961-912	959963-912

Column Brand	Manufacturer	5 µm				
		2.1 x 50 mm	4.6 x 50 mm	4.6 x 100 mm	4.6 x 150 mm	4.6 x 250 mm
Avantor® ACE® Generix® Phenyl Hexyl		GEN-5PHEX-0502	GEN-5PHEX-0546	GEN-5PHEX-1046	GEN-5PHEX-1546	GEN-5PHEX-2546
Acclaim 120 Phenyl-1	Thermo Scientific	-	-	-	-	079697
Ascentis Phenyl	Supelco	581611-U	581615-U	-	581616-U	581617-U
Hypersil Phenyl	Thermo Scientific	-	-	-	30905-154630	30905-254630
Hypersil Phenyl-BDS	Thermo Scientific	-	-	-	28905-154630	28905-254630
Hypersil GOLD Phenyl ¹	Thermo Scientific	-	25905-054630	-	25905-154630	25905-254630
Inertsil Ph	GL Sciences	-	-	-	5020-01324	5020-01328
Inertsil Ph-3 ²	GL Sciences	5020-05112	5020-05142	5020-05144	5020-01920	5020-01921
InertSustain Phenyl ²	GL Sciences	5020-16303	5020-16324	5020-16326	5020-16327	5020-16328
InertSustain Phenylhexyl	GL Sciences	5020-89002	5020-89016	5020-89025	5020-89026	5020-89027
Luna Phenyl-Hexyl ³	Phenomenex	00B-4257-B0	00B-4257-E0	00D-4257-E0	00F-4257-E0	00G-4257-E0
Nucleosil C6H5	Macherey-Nagel	-	-	-	-	720956.46
Nucleodur Phenyl-Hexyl ³	Macherey-Nagel	760583.20	760583.46	760586.46	760588.46	760589.46
Prodigy Phenyl-3 ³	Phenomenex	-	-	00B-4298-E0	00F-4298-E0	00G-4298-E0
Purospher STAR Phenyl ²	Merck	-	-	-	1.51919.0001	1.51918.0001
Roc Phenyl-Hexyl	Restek	-	-	9535515	9535565	9535575
Synchronis Phenyl ⁴	Thermo Scientific	-	-	-	97905-154630	97905-254630
Waters Spherisorb Phenyl	Waters	-	-	PSS830812	PSS830813	PSS830815
YMC-Pack Ph	YMC	PH12S05-05Q1QT	PH12S05-0546WT	PH12S05-1046WT	PH12S05-1546WT	PH12S05-2546WT
Zorbax Eclipse XDB-Phenyl ⁵	Agilent	960967-912	935967-912	-	993967-912	990967-912
Zorbax Eclipse Plus-Phenyl-Hexyl ⁵	Agilent	-	-	959996-912	959993-912	959990-912

¹ 1.9 µm; ² 2 µm; ³ 2.0 mm ID; ⁴ 1.7 µm ID; ⁵ 3.5 µm

Note: Phenyl-Hexyl phases incorporate a C6 alkyl spacer which may result in improved hydrophobicity and increased stability in acidic conditions compared to conventional phenyl bonded phases.

Avantor®, ACE® and Generix® are registered trademarks of Avantor, Inc. or its affiliates. All other trademarks are the property of their respective owners

Avantor® ACE® Generix® SIL Cross Reference Guide

For more information
or advice please
contact our technical
support team at
[chromsupport@
avantorsciences.com](mailto:chromsupport@avantorsciences.com)

ALTERNATIVES TO THE FOLLOWING LEADING AND LEGACY COLUMN BRANDS:

(Manufacturers' original part numbers are shown in the table below. Avantor® ACE® Generix® columns are also available in additional dimensions - please enquire for further details.)

Column Brand	Manufacturer	3 µm		
		2.1 x 50 mm	4.6 x 100 mm	4.6 x 150 mm
Avantor® ACE® Generix® SIL		GEN-3SIL-0502	GEN-3SIL-1046	GEN-3SIL-1546
Ascentis Si	Supelco	581500-U	-	-
Hypersil Silica	Thermo Scientific	30003-052130	-	30003-154630
Hypersil GOLD Silica ¹	Thermo Scientific	-	-	25103-154630
Inertsil SIL	GL Sciences	5020-04212	5020-04244	5020-01701
Luna Silica(2) ²	Phenomenex	00B-4162-B0	00D-4162-E0	00F-4162-E0
Nucleosil SiOH	Macherey-Nagel	-	-	-
Nucleodur SiOH	Macherey-Nagel	-	-	760172.46
Roc Silica	Restek	-	9530315	9530365
Purospher STAR Si	Merck	-	-	-
SunFire Silica	Waters	-	-	-
Supelcosil LC-Si	Supelco	-	-	58981
Synchronis Silica ³	Thermo Scientific	-	-	97003-154630
Ultra Silica	Restek	9100352	9100315	9100365
Waters Spherisorb Silica	Waters	-	PSS832012	PSS832013
YMC-Pack SIL	YMC	SL12S03-05Q1QT	SL12S03-1046WT	SL12S03-1546WT

Column Brand	Manufacturer	5 µm				
		2.1 x 50 mm	4.6 x 50 mm	4.6 x 100 mm	4.6 x 150 mm	4.6 x 250 mm
Avantor® ACE® Generix® SIL		GEN-5SIL-0502	GEN-5SIL-0546	GEN-5SIL-1046	GEN-5SIL-1546	GEN-5SIL-2546
Ascentis Si	Supelco	-	-	-	581512-U	581513-U
Hypersil Silica	Thermo Scientific	-	30005-054630	-	30005-154630	30005-254630
Hypersil GOLD Silica ¹	Thermo Scientific	-	-	-	25105-154630	25105-254630
Inertsil SIL	GL Sciences	5020-04312	5020-04342	5020-04344	5020-01711	5020-01712
Luna Silica(2) ²	Phenomenex	00B-4274-B0	00B-4274-E0	00D-4274-E0	00F-4274-E0	00G-4274-E0
Nucleosil SiOH	Macherey-Nagel	-	-	-	-	720099.46
Nucleodur SiOH	Macherey-Nagel	-	760023.46	-	760012.46	760007.46
Roc Silica	Restek	-	-	9530515	9530565	9530575
Purospher STAR Si	Merck	-	-	-	1.50356.0001	1.50357.0001
SunFire Silica ³	Waters	-	-	-	186003453	186003454
Supelcosil LC-Si	Supelco	-	58236	-	58200-U	58295
Synchronis Silica ³	Thermo Scientific	-	-	97005-104630	97005-154630	97005-254630
Ultra Silica	Restek	9100552	9100555	9100515	9100565	9100575
Waters Spherisorb Silica	Waters	-	PSS830111	PSS083112	PSS830113	PSS830115
YMC-Pack SIL	YMC	SL12S05-05Q1QT	SL12S05-0546WT	SL12S05-1046WT	SL12S05-1546WT	SL12S05-2546WT
Zorbax Rx-Sil	Agilent	-	-	-	883975-901	880975-901

¹1.9 µm; ² 2.0 mm ID; ³ 1.7 µm

Avantor®, ACE® and Generix® are registered trademarks of Avantor, Inc. or its affiliates. All other trademarks are the property of their respective owners

Ordering Information

AVANTOR® ACE® GENERIX® 1.8 µm UHPLC COLUMNS

Columns Dimensions	Particle size	UHPLC/HPLC hardware pressure rated up to 1,000 bar / 15,000 psi		
		C18(2)	C8(2)	Phenyl-Hexyl
2.1 x 50 mm	1.8 µm	GEN-18C182-0502U	GEN-18C82-0502U	GEN-18PHEX-0502U
2.1 x 75 mm	1.8 µm	GEN-18C182-7502U	GEN-18C82-7502U	GEN-18PHEX-7502U
2.1 x 100 mm	1.8 µm	GEN-18C182-1002U	GEN-18C82-1002U	GEN-18PHEX-1002U
3.0 x 50 mm	1.8 µm	GEN-18C182-0503U	GEN-18C82-0503U	GEN-18PHEX-0503U
3.0 x 75 mm	1.8 µm	GEN-18C182-7503U	GEN-18C82-7503U	GEN-18PHEX-7503U
3.0 x 100 mm	1.8 µm	GEN-18C182-1003U	GEN-18C82-1003U	GEN-18PHEX-1003U
Guard cartridges for 2.1-3.0 mm ID UHPLC columns (3/pk)		GEN-C182-GDU	GEN-C82-GDU	GEN-PHEX-GDU
UHPLC guard holder		H0011	H0011	H0011

AVANTOR® ACE® GENERIX® 3 µm HPLC COLUMNS

Columns Dimensions	Particle size	HPLC hardware pressure rated up to 275 bar / 4,000 psi			
		C18(2)	C8(2)	Phenyl-Hexyl	Sil
2.1 x 50 mm	3 µm	GEN-3C182-0502	GEN-3C82-0502	GEN-3PHEX-0502	GEN-3SIL-0502
2.1 x 75 mm	3 µm	GEN-3C182-7502	GEN-3C82-7502	GEN-3PHEX-7502	GEN-3SIL-7502
2.1 x 100 mm	3 µm	GEN-3C182-1002	GEN-3C82-1002	GEN-3PHEX-1002	GEN-3SIL-1002
2.1 x 125 mm	3 µm	GEN-3C182-1202	GEN-3C82-1202	GEN-3PHEX-1202	GEN-3SIL-1202
2.1 x 150 mm	3 µm	GEN-3C182-1502	GEN-3C82-1502	GEN-3PHEX-1502	GEN-3SIL-1502
3.0 x 50 mm	3 µm	GEN-3C182-0503	GEN-3C82-0503	GEN-3PHEX-0503	GEN-3SIL-0503
3.0 x 75 mm	3 µm	GEN-3C182-7503	GEN-3C82-7503	GEN-3PHEX-7503	GEN-3SIL-7503
3.0 x 100 mm	3 µm	GEN-3C182-1003	GEN-3C82-1003	GEN-3PHEX-1003	GEN-3SIL-1003
3.0 x 125 mm	3 µm	GEN-3C182-1203	GEN-3C82-1203	GEN-3PHEX-1203	GEN-3SIL-1203
3.0 x 150 mm	3 µm	GEN-3C182-1503	GEN-3C82-1503	GEN-3PHEX-1503	GEN-3SIL-1503
4.0 x 50 mm	3 µm	GEN-3C182-0504	GEN-3C82-0504	GEN-3PHEX-0504	GEN-3SIL-0504
4.0 x 75 mm	3 µm	GEN-3C182-7504	GEN-3C82-7504	GEN-3PHEX-7504	GEN-3SIL-7504
4.0 x 100 mm	3 µm	GEN-3C182-1004	GEN-3C82-1004	GEN-3PHEX-1004	GEN-3SIL-1004
4.0 x 125 mm	3 µm	GEN-3C182-1204	GEN-3C82-1204	GEN-3PHEX-1204	GEN-3SIL-1204
4.0 x 150 mm	3 µm	GEN-3C182-1504	GEN-3C82-1504	GEN-3PHEX-1504	GEN-3SIL-1504
4.6 x 50 mm	3 µm	GEN-3C182-0546	GEN-3C82-0546	GEN-3PHEX-0546	GEN-3SIL-0546
4.6 x 75 mm	3 µm	GEN-3C182-7546	GEN-3C82-7546	GEN-3PHEX-7546	GEN-3SIL-7546
4.6 x 100 mm	3 µm	GEN-3C182-1046	GEN-3C82-1046	GEN-3PHEX-1046	GEN-3SIL-1046
4.6 x 125 mm	3 µm	GEN-3C182-1246	GEN-3C82-1246	GEN-3PHEX-1246	GEN-3SIL-1246
4.6 x 150 mm	3 µm	GEN-3C182-1546	GEN-3C82-1546	GEN-3PHEX-1546	GEN-3SIL-1546
Guard cartridges for 2.1 mm ID HPLC columns (5/pk)		GEN-C182-GD2	GEN-C82-GD2	GEN-PHEX-GD2	GEN-SIL-GD2
Guard cartridges for 3.0-4.6 mm ID HPLC columns (5/pk)		GEN-C182-GD	GEN-C82-GD	GEN-PHEX-GD	GEN-SIL-GD
HPLC guard holder		H0010	H0010	H0010	H0010

AVANTOR® ACE® GENERIX® 5 µm HPLC COLUMNS

Columns Dimensions	Particle size	HPLC hardware pressure rated up to 275 bar / 4,000 psi			
		C18(2)	C8(2)	Phenyl-Hexyl	Sil
2.1 x 50 mm	5 µm	GEN-5C182-0502	GEN-5C82-0502	GEN-5PHEX-0502	GEN-5SIL-0502
2.1 x 75 mm	5 µm	GEN-5C182-7502	GEN-5C82-7502	GEN-5PHEX-7502	GEN-5SIL-7502
2.1 x 100 mm	5 µm	GEN-5C182-1002	GEN-5C82-1002	GEN-5PHEX-1002	GEN-5SIL-1002
2.1 x 125 mm	5 µm	GEN-5C182-1202	GEN-5C82-1202	GEN-5PHEX-1202	GEN-5SIL-1202
2.1 x 150 mm	5 µm	GEN-5C182-1502	GEN-5C82-1502	GEN-5PHEX-1502	GEN-5SIL-1502
2.1 x 250 mm	5 µm	GEN-5C182-2502	GEN-5C82-2502	GEN-5PHEX-2502	GEN-5SIL-2502
3.0 x 50 mm	5 µm	GEN-5C182-0503	GEN-5C82-0503	GEN-5PHEX-0503	GEN-5SIL-0503
3.0 x 75 mm	5 µm	GEN-5C182-7503	GEN-5C82-7503	GEN-5PHEX-7503	GEN-5SIL-7503
3.0 x 100 mm	5 µm	GEN-5C182-1003	GEN-5C82-1003	GEN-5PHEX-1003	GEN-5SIL-1003
3.0 x 125 mm	5 µm	GEN-5C182-1203	GEN-5C82-1203	GEN-5PHEX-1203	GEN-5SIL-1203
3.0 x 150 mm	5 µm	GEN-5C182-1503	GEN-5C82-1503	GEN-5PHEX-1503	GEN-5SIL-1503
3.0 x 250 mm	5 µm	GEN-5C182-2503	GEN-5C82-2503	GEN-5PHEX-2503	GEN-5SIL-2503
4.0 x 50 mm	5 µm	GEN-5C182-0504	GEN-5C82-0504	GEN-5PHEX-0504	GEN-5SIL-0504
4.0 x 75 mm	5 µm	GEN-5C182-7504	GEN-5C82-7504	GEN-5PHEX-7504	GEN-5SIL-7504
4.0 x 100 mm	5 µm	GEN-5C182-1004	GEN-5C82-1004	GEN-5PHEX-1004	GEN-5SIL-1004
4.0 x 125 mm	5 µm	GEN-5C182-1204	GEN-5C82-1204	GEN-5PHEX-1204	GEN-5SIL-1204
4.0 x 150 mm	5 µm	GEN-5C182-1504	GEN-5C82-1504	GEN-5PHEX-1504	GEN-5SIL-1504
4.0 x 250 mm	5 µm	GEN-5C182-2504	GEN-5C82-2504	GEN-5PHEX-2504	GEN-5SIL-2504
4.6 x 50 mm	5 µm	GEN-5C182-0546	GEN-5C82-0546	GEN-5PHEX-0546	GEN-5SIL-0546
4.6 x 75 mm	5 µm	GEN-5C182-7546	GEN-5C82-7546	GEN-5PHEX-7546	GEN-5SIL-7546
4.6 x 100 mm	5 µm	GEN-5C182-1046	GEN-5C82-1046	GEN-5PHEX-1046	GEN-5SIL-1046
4.6 x 125 mm	5 µm	GEN-5C182-1246	GEN-5C82-1246	GEN-5PHEX-1246	GEN-5SIL-1246
4.6 x 150 mm	5 µm	GEN-5C182-1546	GEN-5C82-1546	GEN-5PHEX-1546	GEN-5SIL-1546
4.6 x 250 mm	5 µm	GEN-5C182-2546	GEN-5C82-2546	GEN-5PHEX-2546	GEN-5SIL-2546
Guard cartridges for 2.1 mm ID HPLC columns (5/pk)		GEN-C182-GD2	GEN-C82-GD2	GEN-PHEX-GD2	GEN-SIL-GD2
Guard cartridges for 3.0-4.6 mm ID HPLC columns (5/pk)		GEN-C182-GD	GEN-C82-GD	GEN-PHEX-GD	GEN-SIL-GD
HPLC guard holder		H0010	H0010	H0010	H0010

AVANTOR® ACE® GENERIX® HPLC COLUMN MULTIPACKS

Three column multipacks provide further cost savings

Column Dimensions	Particle size	HPLC hardware pressure rated up to 275 bar / 4,000 psi			
		C18(2)	C8(2)	Phenyl-Hexyl	Sil
2.1 x 50 mm (3/pk)	3 µm	GEN-3C182-0502/3PK	GEN-3C82-0502/3PK	GEN-3PHEX-0502/3PK	GEN-3SIL-0502/3PK
4.6 x 100 mm (3/pk)	3 µm	GEN-3C182-1046/3PK	GEN-3C82-1046/3PK	GEN-3PHEX-1046/3PK	GEN-3SIL-1046/3PK
4.6 x 150 mm (3/pk)	3 µm	GEN-3C182-1546/3PK	GEN-3C82-1546/3PK	GEN-3PHEX-1546/3PK	GEN-3SIL-1546/3PK
2.1 x 50 mm (3/pk)	5 µm	GEN-5C182-0502/3PK	GEN-5C82-0502/3PK	GEN-5PHEX-0502/3PK	GEN-5SIL-0502/3PK
4.6 x 150 mm (3/pk)	5 µm	GEN-5C182-1546/3PK	GEN-5C82-1546/3PK	GEN-5PHEX-1546/3PK	GEN-5SIL-1546/3PK
4.6 x 250 mm (3/pk)	5 µm	GEN-5C182-2546/3PK	GEN-5C82-2546/3PK	GEN-5PHEX-2546/3PK	GEN-5SIL-2546/3PK

Setting science in motion to create a better world

Avantor® is a leading global provider of mission critical products and services to customers in the biopharma, healthcare, education & government, and advanced technologies & applied materials industries. We operate in more than 30 countries and deliver an extensive portfolio of products and services.